



HEALTHCARE SEATING

6200 SERIES

PHLEBOTOMY CHAIR

The Lumex 6200 Series bariatric blood draw/phlebotomy chair is designed with both the patient and phlebotomist in mind. The sensible and affordable design offers a small footprint and modern style while keeping patient comfort and well-being the highest priority.

Features:

- Fully upholstered padded seat, back, and arm rests
- Shown with optional wall saver wireframe
- Powder coat frame finish - Foghorn
- Available in 24 industry standard colors to complement most decors
- Individual height adjustable arm rests

[6200 Product Info](#) [Fabric Choices](#)



Single Flip Arm
Dual Flip Arm

6200 Options:

6200-xx	■
6200DA-xx	■

★ options available at additional cost

6200 Accessories:

Wallsaver Wireframe [§]	item no. 620-WALL SAVER
----------------------------------	----------------------------

accessories sold separately at additional cost



Specifications:

Series	Overall Width		Overall Depth		Overall Height		Seat Width		Seat Depth		Seat Height		Armrest Height		Back Width	Back Height	Shipping Weight [§]	Weight Capacity ¹						
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	low	high				lb	kg [†]	lb	kg [†]			
6200	34.00	86.36	29.00	73.66	37.50	95.25	24.50	62.23	18.50	46.99	22.00	64.77	27.00	62.23	36.00	66.04	25.50	64.77	17.00	43.18	55	24.95	800	362.87

★ options available at additional cost †-xx in the product number is the placeholder for fabric choice § images show this accessory in use - sold separately at additional cost
 † standard rounding method used ¹ maximum weight capacity - evenly distributed - refer to user manual § approximate measurement - options alter spec



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF1800043 Rev1/25