

BARIATRIC/HEAVY DUTY WHEELCHAIRS

TRAVELER® HD WHEELCHAIR

The Traveler® HD is manufactured with the same commitment to quality E&J is known for.

Designed for optimal patient transport as well as the hospital/clinical use. Built with a heavy duty frame to support up to 500 lb users. Reinforced at all critical load-bearing points and a long lasting embossed back and seat upholstery with a double inner liner for durability. With its multiple seat widths ranging from 20" to 24" and 18" seat depth (1" longer than standard), this heavy duty wheelchair enables institutions and caregivers to accommodate a large variety of patients.

Features

- Available in 20", 22" and 24" seat widths, with 18" seat depth (1" longer than standard)
- Carbon steel reinforced frame include double cross braces, heavy-duty 5/8" diameter axles
- Embossed black leatherette upholstery with double inner lining for durability with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Two-position axle and caster fork provide 2" seat height adjustment for a hemi position
- Tool-free adjustable front riggings with durable aluminum black footplates
- Rear 24" x 1" black plastic wheels with black handrims and non-marking polyurethane tires
- Maximum Weight Capacity 500 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007

Warranty:

Five (5) year limited warranty



with Elevating Legrests



Anti-fold bar



Wheel-lock



desk arms (shown) and
swingaway footrest (shown)
— other configurations available



Color variations may occur due to monitor or printer settings.

