

REPLACEMENT PARTS

FR588W POWER BARIATRIC
RECLINER - WIDE

Item	#	UOM	Description
5792A001G	1	1 ea	Table, Side
55792A002	2	4 ea /set	Hardware, Clevis Pin w Washer, Side Table
58806U(COLOR CODE)FR	3	1 ea	Headrest w/ Flap
58805U(COLOR CODE)FR	4	1 ea	Armrest, Upholstered, Left
58804U(COLOR CODE)FR	5	1 ea	Armrest, Upholstered, Right
58803U(COLOR CODE)FR	6	1 ea	Back, Upholstered
58802U(COLOR CODE)FR	7	1 ea	Footrest, Upholstered
58801U(COLOR CODE)FR	8	1 ea	Seat, Upholstered
5880A122FR	10	1 pr	Side Panels
AP020004-8	11	4 ea /set	Hardware, Push Pin & Spring Side Panels
588W016A	12	1 ea	Footrest
FR588W016W	14	1 ea	Armrest Support, Left (not shown)
FR588W011W	15	1 ea	Armrest Support, Right (not shown)
588W016W	16	1 ea	Armrest Support, Left (not shown)
588W011W	17	1 ea	Armrest Support, Welded, Right (not shown)



Item	#	UOM	Description
WB170088	18	1 ea	Caster, Rear w/ Lock
WC170088	19	1 ea	Caster, Front
588WZA001	20	1 ea	Trigger Release Assembly
EI-000002	21	1 ea	Hardware, Gas Spring
LRMP0001	22	1 ea	Paint Touch-Up (Lumex Gray)
588BKACTIONATOR	24	1 ea	Actuator, Back / Seat Section (not shown)
588CNTRLBX	25	1 ea	Actuator, Control Box
588CNTRLBXCABLE	25	1 set	Cables for Control Box (not shown)
588CNTRBXPWRCBL	25	1 set	Power Cable Set (not shown)
588HANDSET	26	1 ea	Actuator, Handset (New)
HB-72600001	27	1 ea	Actuator, Handset- Old Style (not shown)
588LGACTUATOR	28	1 ea	Actuator, Legrest Section (not shown)
A011F	29	1 ea	Battery Back, Actuator (not shown)

TO ORDER: Tel: 1.770.368.4700 • www.grahamfield.com

Returns: Parts may not be returned without prior authorization from GF Health Products, Inc. Custom parts cannot be returned for credit. Please refer to GF Health Products, Inc. [Terms and Conditions](#) for further information.



GF Health Products, Inc. is an
Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a trademark of GF Health Products, Inc.
GF Health Products, Inc. is an ISO 13485 Certified Company.

**graham
field**