

REPLACEMENT PARTS

FR601 PIVOT ARM CLINICAL CARE



Item #	UOM	Description	Item #	UOM	Description
5792A001G	1 ea	Table, Side	CSTRTRGRN01	1 ea	Caster, Directional, Green Tab Tente
5792A002	4 ea /set	Hardware, Clevis Pin w Washer, Side Table	CSTRTRBLK01	1 ea	Caster, Swivel, Locking Black Tab Tente
FR601G030M	1 ea	Trendelenburg Release Bar	601G008W	1 ea	Pedal, Swing Arm Release
PV860122	1 ea	Rubber Grip, Trendelenburg Bar	5772A216G	1 ea	Push Handle w Grip
601U09J(color code)FR	1 ea	Pivot Arm, Upholstered w Arm Cap, Right	9205A	1 ea	Push Handle, Grip Only
601P09J(color code)FR	1 ea	Pivot Arm, Upholstered, Right	PP46000	1 ea	Hardware, Back Attachment Thumb Screws
601U10J(color code)FR	1 ea	Pivot Arm, Upholstered w Arm Cap, Left	5660A200	1 ea	Mechanism, Recliner
601P10J(color code)FR	1 ea	Pivot Arm, Upholstered, Left	LRHM0114	1 ea	Pendant Holster, Heat & Massage
601SPFTPLATE	1 ea	Footplate, Self-Storing	HM00G3CORD	1 ea	Powercord Supply, Heat & Massage Gen 3
60102J(color code)FR	1 ea	Back, Upholstered	ECMCGF2	1 ea	Control Box, Heat & Massage Gen 3
60102JH(color code)FR	1 ea	Back, Upholstered, w Heat & Massage	HM00G3PEND	1 ea	Pendant, Heat & Massage Gen 3
60102JHO(color code)FR	1 ea	Back, Upholstered, Heat Only	HM00G3PENDHO	1 ea	Pendant, Heat Only Gen 3
60101J(color code)FR	1 ea	Seat, Upholstered	LRMP0001	1 ea	Paint, Touch-Up, Gray (not shown)
60101JH(color code)FR	1 ea	Seat, Upholstered, Heat & Massage / Heat Only			
60104J(color code)FR	1 ea	Ottoman, Upholstered			
60106J(color code)FR	1 ea	Footrest, Upholstered			



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company



Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2023 GF Health Products, Inc. All Rights Reserved.
Lumex is a trademark of GF Health Products, Inc.
GF Health Products, Inc. is an ISO 13485 Certified Company.

graham
field