




REPLACEMENT PARTS

FR587W DELUXE CLINICAL CARE - WIDE



Heat and Massage Feature —
Control Module and Transformer Unit
Showing All Connections Correctly Made

Item #	UOM	Description
5792A001G	1 ea	Table, Side
5792A002	4 ea /set	Hardware, Clevis Pin w/ Washer, Side Table
58737J(color code)FR	1 ea	Ottoman, Upholstered
58735J(color code)FR	1 ea	Armrest, Upholstered
58733J(color code)FR	1 ea	Footrest, Upholstered
58734J(color code)FR	1 ea	Legrest, Upholstered
58731J(color code)FR	1 ea	Seat, Upholstered
58731JH(color code)FR	1 ea	Seat, Upholstered, Heat & Massage / Heat Only
56605J(color code)FR	1 ea	Side Panels, Upholstered (not shown)
58732J(color code)FR	1 ea	Back, Upholstered
58732JH(color code)FR	1 ea	Back, Upholstered, Heat & Massage
58732JHO(color code)FR	1 ea	Back, Upholstered, Heat Only
AP020004-8	8 ea / set	Hardware, Push Pin and Spring, Side Panels
CSTRTRGN01	1 ea	Caster, Directional Green Tab Tente

Item #	UOM	Description
CSTRBLK01	1 ea	Caster, Swivel, Locking Black Tab Tente
5873C400A	1 ea	Push Handle w Grip
9205	1 ea	Grip Only, Push Handle
5610F200	1 ea	Mechanism, Recliner
5770A073	1 ea	Hardware, Gas Spring Bracket
EI000001	1 ea	Hardware, Gas Spring
5792A003	10 ea /set	Hardware, Clevis Pin w Cotter Clip
LRHM0114	1 ea	Pendant Holster, Heat & Massage
HM00G3CORD	1 ea	Powercord Supply, Heat & Massage Gen 3
ECMCGF2	1 ea	Control Box, Heat & Massage Gen 3
HM00G3PEND	1 ea	Pendant, Heat & Massage Gen 3
HM00G3PENDHO	1 ea	Pendant, Heat Only Gen 3
5772A001	1 ea	Hardware, Cable Assemble (not shown)
LRMP0001	1 ea	Paint, Touch-Up, Gray (not shown)



GF Health Products, Inc. is an
Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

•1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a trademark of GF Health Products, Inc.
GF Health Products, Inc. is an ISO 13485 Certified Company.

graham
field