

WARMERS

TABLE TOP WARMERS

The tabletop warming cabinets come in two sizes, standard and compact. It is suggested that these units sit on an existing cabinet or on a Blickman mobile stand. One adjustable perforated stainless steel shelf comes standard with the standard unit and is available as an option for the compact unit. The door can come in either glass or stainless steel.

Features

- Two sizes: standard and compact
- Standard Electrical
 - Amps: 7.5
 - Total Wattage: 750
- Compact Electrical
 - Amps: 3.5
 - Total Wattage: 350
- Size appropriate for the tabletop location
- Standard version includes a lock to secure contents
- One perforated stainless steel shelf comes standard
- Ability to set to solutions-only warming cabinet by programming to heat only between 90°F and 110°F (32.2°C - 43.3°C)
- Heat variances between 90°F and 160°F (32.2°C - 71.1°C)
- The temperature remains within 1 degree Fahrenheit of set point with virtually no leakage into adjoining chambers or outside area
- The fan is located at the top, allowing for better air circulation, and frees up the bottom shelf for more storage
- Flat Top



		overall Width		overall Height		overall Depth		UOM
		inch	cm	inch	cm	inch	cm	
GF-14B7925243	Table Top Warmer - Glass Door	30	76.2	24.5	62.2	26.625	67.6	1 ea
GF-14B7925200	Table Top Warmer - Solid Door	30	76.2	24.5	62.2	26.625	67.6	1 ea
GF-14B7927243	Table Top Warmer - Compact, Glass Door	24	61	24.5	62.2	20.625	52.4	1 ea
GF-14B7927200	Table Top Warmer - Compact, Solid Door	24	61	24.5	62.2	20.625	52.4	1 ea



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2025, GF Health Products, Inc. All Rights Reserved. Graham-Field is a trademark of GF Health Products, Inc.