

WARMERS

UNDERCOUNTER WARMERS

This warming cabinet is suited for use under a counter-top or in limited spaces. Its door comes in either glass or stainless steel. Each unit comes standard with two adjustable perforated stainless steel shelves. Each unit comes standard with two adjustable perforated stainless steel shelves.

Features

- Size appropriate for a location under a counter
- One glass and stainless steel door, lock to secure contents
- Two perforated stainless steel shelves comes standard
- Electrical
 - Amps: 7.5
 - Total Wattage: 750
- Heat variances between 90°F and 160°F (32.2°C - 71.1°C)
- The temperature remains within 1 degree Fahrenheit of set point with virtually no leakage into adjoining chambers or outside area
- The fan is located at the top, allowing for better air circulation, and frees up the bottom shelf for more storage
- Flat top



		overall Width		overall Height		overall Depth		UOM
		inch	cm	inch	cm	inch	cm	
GF-14B7922243	Undercounter Warmer - Glass Door	30	76.2	35.5	90.2	26.625	67.6	1 ea
GF-14B7922200	Undercounter Warmer - Solid Door	30	76.2	35.5	90.2	26.625	67.6	1 ea



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2025, GF Health Products, Inc. All Rights Reserved. Graham-Field is a trademark of GF Health Products, Inc.