

# WARMERS

## DUAL CHAMBER FREESTANDING WARMERS

This freestanding cabinet includes two chambers with doors that come in either glass or stainless steel and lock to secure contents. Each unit comes standard with three adjustable perforated stainless steel shelves for flexibility and expanded use, and integrated casters for ease of mobility.

### Features

- Freestanding
- Dual-chamber model
- Two glass or stainless steel doors, lock to secure contents
- Self-closing, magnetic gasket, adjustable tension doors afford convenient hands-free opening and closing
- Three perforated stainless steel shelves come standard and are adjustable for flexibility and expanded use
- Electrical
  - Amps: 15
  - Total Wattage: 1500
- Heat variances between 90°F and 160°F (32.2°C - 71.1°C)
- The temperature remains within 1 degree Fahrenheit of set point with virtually no leakage into adjoining chambers or outside area
- The fan is located at the top, allowing for better air circulation, and frees up the bottom shelf for more storage
- Flat top
- Integrated, hidden 3" (7.6 cm) ball-bearing polyolefin casters with a zinc finish for ease of mobility



		overall Width		overall Height		overall Depth		UOM
		inch	cm	inch	cm	inch	cm	
GF-14B792423	Dual Chamber Freestanding Warming Cabinet - Glass Door	30	76.2	74.5	189.2	26.625	67.6	1 ea
GF-14B7924200	Dual Chamber Freestanding Warming Cabinet - Solid Door	30	76.2	74.5	189.2	26.625	67.6	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2025, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.