

WHEELCHAIR ACCESSORIES

WHEELCHAIR POSITIONING TRAY

Multipurpose tray is ideal for eating, reading, and writing while seated in wheelchair or resting in bed

Features

- Ergonomic design with contoured cutout and comfortable beanbag style base
- Easy to clean composite top measures 15.75" x 13" x 2" including 2.5" diameter (40 cm x 33 cm x 5 cm including 6.5 cm diameter) engraved circumferences for sitting cups or bottles
- Washable nylon cover includes interior 6" (15.24 cm) zippered small pocket on top and side pocket
- Convenient carry handle attached
- Composite tray can be cleaned with soap and water or diluted bleach
- Cover can be removed and hand washed with soap and water, re-attached to fastener tabs
- Total weight: 1.3 lb (600g)
- Weight capacity: 10 lb (4.53kg) evenly distributed
- This product is not made with natural rubber latex

Warranty

One (1) year limited warranty



	overall Width		overall Depth		overall height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	mm	inch	mm	inch	mm	lb	g	lb	kg	
EJTRAY18	15.75	40	13	33	2	5.1	1.3	600	10	4.53	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham
field

GF1900048RevC25