

WHEELCHAIR LEGRESTS/FOOTRESTS

SWINGAWAY FOOTRESTS WITH COMPOSITE FOOTPLATE

Features

- "Tool Free" length adjustment
- Light weight composite footplate
- 1.5" (3.8 cm) industry standard spacing for attachment
- Long lasting and durable hammertone finish
- Includes heel loop to ensure foot position
- Latex free
- Standard equipment on E&J Advantage, Traveler SE, Traveler L3 and Traveler L4
- Meets Medicare code K0052



		UOM
90763030	Composite Footplate	1 pr



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2014, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.