



WHEELCHAIR ACCESSORIES

PNEUMATIC TIRE
AND WHEEL KIT

Features

- 24" x 1 3/8" gray non-marking nylon tire
- Comes standard with 7/16" bearings
- Kit also includes 12 mm bearings
- Black composite wheel with protective hub cap included
- Black composite handrim
- Air inflate to 50 lb
- Intended for use on chairs up to 250 lb weight capacity - EVENLY DISTRIBUTED



		tire Diameter		tire Width		tire Pressure		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	psi	bar	lb	kg	
90763056	Gray	24	61	1.375	3.5	50	3.45	250	113.4	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2014, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham
field

GF1400026RevB35