

RESPIRATORY

JB0112-067 NEBULIZER COMPRESSOR NEB-U-TYKE® SPEEDSTER

This interactive, child-friendly design provides an ideal solution for effective aerosol treatment for children. Designed with an attractive car shape, this will keep children relaxed and entertained during aerosol treatments. The Neb-u-Tyke® Speedster is the latest addition to the Graham-Field line of quality nebulizer compressors that offer innovation, performance, and value for both patients and providers.

Features:

- Light weight, reliable, piston-driven compressor
- Colorful, child-friendly design complete with stickers to personalize during aerosol treatment
- Convenient and compact design for easy transport and storage
- Comes complete with Nebulizer Accessory Kit
- Designed to meet the requirements of HCPCS Code: E0570
- Five year limited warranty



Package Contains:

- 1 Neb-u-Tyke® Speedster Nebulizer Compressor
- 1 User Manual
- 1 "T" Adaptor
- *1 Nebulizer Accessory Kit *complete with*
 - 1 Reusable nebulizer (*single patient use only*)
 - 1 Pediatric mask
 - 1 Adult Mask
 - 1 Mouthpiece
 - 1 PVC tubing
 - 3 Replacement filters



Options:

★ Item No. UOM Description

Item No.	UOM	Description
JB0112-067	1 ea	Neb-u-Tyke® Speedster Nebulizer
JB77200	1 ea	Neb-u-Tyke® Speedster with Garage
JB76220	1 ea	Garage Bag Only

†Specifications:

Overall Width		Overall Length		Overall Height		Operating Range				Storage Range				Power		Product Weight		Sound Level	Max Pressure	Max Flow	Particle Size	Capacity				
low	high	low	high	low	high	Temperature		Humidity		Temperature		Humidity		VAC	Hz	lb	kg	dBA	compressor	nebulizer						
inch	cm	inch	cm	inch	cm	°F	°C	°F	°C	RH	RH	°F	°C	°F	°C	RH	RH		psi	lpm	µm	ml				
6.9	17.52	10.2	25.9	3.7	9.39	41	4.4	104	40	30%	95%	14	-10	158	70	10%	95%	115	60	3.8	1.72	<55	30	≈12	≤5	5

*Boiling nebulizer components is not recommended.

★FDA (USA) restricts the sale of this device to or by the order of a physician

†Specifications are subject to change without notice.



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Neb-u-Tyke is a registered trademark of GF Health Products, Inc.