

RESPIRATORY

JB0112-016 NEBULIZER COMPRESSOR VACUTEC™ 800 EV2

The Vacutec™ 800 EV2 Aspirator is a professional medical suction device that produces a maximum vacuum of 560 mmHg. The unit is equipped with an 800 cc collection bottle and has been designed for convenient and reliable use for most home settings.

Features:

- Designed for most homecare suction applications
- Low level noise at less than 55 dBA
- Designed for single patient, everyday use
- Easy-to-read gauge and easy-to-turn flow regulator knob
- Suction cup feet for stable operation
- Complete with disposable 800cc collection bottle
- Automatic shut-off
- Complies with EN60601 standards
- Designed to meet the requirements of HCPCS Code: E0570
- Two year limited warranty

Package Contains:

- 1 Vacutec™ 800 EV2 Aspirator
- 1 User Manual
- 1 Collection bottle
- 1 Connection tubing (from aspirator to collection bottle)
- 1 Pack of suction tubing (one-time patient use)
- 1 Replacement HEPA filter

Options:

★ Item No.	UOM	Description
JB0112-016	1 ea	Vacutec™ 800 EV2 Aspirator
JB0112-016B	6/cs	Collection Bottle
JB0112-016F	3 ea/pk	Replacement HEPA Filter
JB0112-016S	1 ea	Suction Tubing - 6 ft long
JB0112-016C	1 ea	Connection Tubing



†Specifications:

†Specifications:																														
Overall Width		Overall Length		Overall Height		Operating Range						Storage Range						Power		Product Weight		Sound Level	Consumption		Flow	compressor		Vacuum Range		Power Consumption
						Temperature		Humidity		Temperature		Humidity		watts	lpm															
						low	high	low	high	low	high	low	high	bacteria filter	mm Hg															
inch	cm	inch	cm	inch	cm	°F	°C	°F	°C	RH	RH	°F	°C	°F	°C	RH	RH	VAC	Hz	lb	kg	dBA	watts	lpm	HEPA	0-560	watts			
4.9	12.44	15.3	38.8	10	25.4	32	5	104	40	30%	95%	14	-10	158	70	10%	95%	115	60	11.6	5.26	<55	60	0-40	HEPA	0-560	60			

★FDA (USA) restricts the sale of this device to or by the order of a physician

†Specifications are subject to change without notice.



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com



© 2013, GF Health Products, Inc. All Rights Reserved.
Neb-u-Tyke is a registered trademark of GF Health Products, Inc.