

RESPIRATORY

NEBULIZER COMPRESSOR NEB-U-LITE EV2

The Neb-u-Lite EV2 is a great addition to the Graham-Field line of quality nebulizer compressors that offer innovation, performance, and value for both patients and providers.

Features:

- Small footprint at a competitive price point for providers
- Light weight Compressor 2.7lb (1.2 kg)
- Built-in nebulizer cup holder for hands-free use
- Reliable, piston-driven compressor is efficient and delivers consistent particle size of <5 microns
- Powerful piston pump is efficient, dependable and quiet
- Complete with disposable Nebulizer Accessory Kit
- Designed to meet the requirements of HCPCS Code: E0570
- Five year limited warranty

Package Contains:

- 1 Neb-U-Lite EV2 Compressor
- *1 Nebulizer Accessory Kit
 - 1 disposable nebulizer cup (reusable for single patient use only)
 - 1 Mouthpiece
 - 1 PVC tubing
 - 5 Spare filters
 - 1 Adult and 1 Pediatric Mask

Options:

★ Item No. UOM Description

JB0112-090	10 ea / cs	Neb-U-Lite EV2
------------	------------	----------------



†Specifications:

Overall Width		Overall Length		Overall Height		Operating Range				Storage Range				Power		Product Weight		Sound Level	Max Pressure	Max Flow	Particle Size Capacity							
						Temperature		Humidity		Temperature		Humidity																
						low	high	low	high	low	high	low	high															
inch	cm	inch	cm	inch	cm	°F	°C	°F	°C	RH	RH	°F	°C	°F	°C	RH	RH	VAC	Hz	lb	kg	dBA	compressor	psi	lpm	nebulizer	µm	ml
5.8	14.73	5.5	13.97	3.3	8.38	50	10	104	40	30%	95%	-40	-40	158	70	10%	95%	120	60	2.7	1.2	<60	42	≈ 8	≤5	8		

*Boiling nebulizer components is not recommended

★FDA (USA) restricts the sale of this device to or by the order of a physician

†Specifications are subject to change without notice.



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

•1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved.
Neb-U-Lite is a registered trademark of GF Health Products, Inc.