



MOBILITY ACCESSORIES

FIXED WHEELS

Features

- Converts standard walker into wheeled walker
- Provides forward mobility with no lateral movement
- The large 5" wheels are ideal for navigating difficult terrain such as pavement, grass or carpet
- Now available with eight adjustment holes to accommodate a broader range of user heights
- Fixed wheels attach to the front of walkers only and are recommended for use with Lumex glide brakes
- Comes complete with glide covers
- Designed to fit 1" standard walker tubing
- Designed to meet the requirements of HCPCS code: E0155
- Limited lifetime warranty



*These wheels are recommended for use on Lumex Walkers.
(We cannot guarantee the performance of these accessories on other industry products.)*

Specifications

Item	Walker Wheels	UOM	Wheel Diameter		Capacity	
			inch	cm [†]	lb	kg [†]
603630A	Small Fixed	pair	3	7.62	300	136.08
603650A	Large Fixed	pair	5	12.70	300	136.08
603750A	Large Fixed with Tall Leg Extensions	pair	5	12.70	300	136.08
603950A	Bariatric Dual Large Fixed	pair	5	12.70	600	272.16

[†] standard rounding method used [†] maximum weight capacity - evenly distributed - refer to user manual



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

*1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.