



SUPPORT SURFACES



CONVERTIBLE MATTRESS

LCM CONVERTIBLE MATTRESS SERIES

The Lumex Convertible Series Mattress provides support to the patient for pressure redistribution. The Lumex mattresses are specifically developed with patient comfort and pressure injury management in mind. Provides pressure reduction and redistribution

Features

- Provides outstanding pressure displacement and redistribution to assist in the prevention and treatment of stage I through IV pressure ulcers

Stretch Cover features

- Air-permeable (breathable)
- Reduces friction and patient shear
- Fluid resistant
- Concealed full zipper
- Easy to clean and maintain

Mattress Construction

- Closed System (Static Air) can quickly convert to an Alternating Pressure Mattress
- Eight non-circular cylinders prevent "rolling" or "gapping" and are fully enclosed in a fiber mesh sleeve and jacket
- 2.5" Hybrid Visco-elastic top layer does not raise body heat as with standard visco-elastic foam
- Air control unit is simple, quiet and set at the optimum 10 minute cycle to eliminate guesswork for the caregiver
- Available with 2" raised side bolsters

Use of the Lumex Convertible will eliminate the need for:

- Costly rental and storage of mattresses that are only used on an as needed basis
- Moving and transporting mattresses from room to room
- Overlays that are low in profile and do not provide adequate support or therapeutic benefit
- Designed to meet the requirements of HCPCS Code E0184



Specifications

Item	UOM	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Shipping ⁶ Weight		Weight Capacity ⁴ EVENLY DISTRIBUTED		Raised Side Bolsters	Defined Perimeter	Soft Heel Zone	Soft Heel Slope	Pressure Reduction	Zippered Cover	Meets 16CFR1633	Made in USA	Mattress in months	Cover in months	
			inch	cm	inch	cm	inch	cm	lb	kg	lb	kg											warranty
LCM3580	1/ea	E0184	36 [†]	91.44	80	203.20	7	17.78	30	13.61	500	226.80	—	■	—	■	■	■	■	■	■	60	12
LCM3580RB	1/ea	E0184	36 [†]	91.44	80	203.20	7	17.78	30	13.61	500	226.80	■	■	—	■	■	■	■	■	■	60	12
LCM4280	1/ea	E0184	42	106.68	80	203.20	7	17.78	30	13.61	500	226.80	—	■	—	■	■	■	■	■	■	60	12
LCM4280RB	1/ea	E0184	42	106.68	80	203.20	7	17.78	30	13.61	500	226.80	■	■	—	■	■	■	■	■	■	60	12
LCM4280B	1/ea	E0184	42	106.68	80	203.20	7	17.78	30	13.61	750	340.19	—	■	—	■	■	■	■	■	■	60	12
LCM4880	1/ea	E0184	48	121.92	80	203.20	7	17.78	35	15.88	500	226.80	—	■	—	■	■	■	■	■	■	60	12
LCM4880RB	1/ea	E0184	48	121.92	80	203.20	7	17.78	35	15.88	500	226.80	■	■	—	■	■	■	■	■	■	60	12
LCM4880B	1/ea	E0184	48	121.92	80	203.20	7	17.78	35	15.88	750	340.19	—	■	—	■	■	■	■	■	■	60	12
LCMP100	1/ea	E0184	5.51	14.00	11.20	28.45	2.76	7.01	9.5	4.31	500	226.80	—	—	—	—	■	—	—	■	pump	12	

■ Standard
 ■ Optional
 ■ Available
 ■ Not Available
 ■ Accessory
 — Does Not Apply
 ■ Active Hyperlink
[†] ± 1 inch tolerance
¹ maximum weight capacity - evenly distributed
⁸ approximate - as size changes so does weight

GF Health Products, Inc. is an
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Product Hyperlinks Enabled
+1.770.368.4700

Information contained herein is subject to change without notice.
 The most current and complete product information can be found on our website.

www.grahamfield.com



© 2016, GF Health Products, Inc. All Rights Reserved.
 Lumex is a registered trademark of GF Health Products, Inc.



GF-1600115 Rev C25