

WHEELCHAIR ACCESSORIES

WHEELCHAIR ONE-HOOK I.V. POLE

**Features**

- Attaches to Advantage<sup>®</sup>, Traveler<sup>®</sup> SE, L3 and old Chromed Traveler Recliner only
- Sure-grip clutch nut provides full range of height positions (maximum height from base rail of chair is 75 inches)
- Designed to meet the requirements of HCPCS code: E0776

**Warranty**

One (1) year limited warranty



		UOM
GF1844W	Chrome	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2025, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.