

BATH ACCESSORIES

SURE-SAFE<sup>®</sup>  
BATHTUB SAFETY RAIL

**Features**

- Heavy-gauge, welded steel construction
- Clamp holds firmly to tub wall without marring tub surface
- Provides two gripping surfaces to facilitate both entry and exit from tub
- Color: White
- Fits tubs with walls 2 3/4" - 6 1/2" thick
- Not recommended for use with fiberglass tubs
- Tool-free installation
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



	Height		UOM
	inch	cm	
2032A	20	50.8	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.