



EXTENDED CARE BEDS



MATRIX SERIES



The Matrix bed features sealed ball bearings in the joints for smooth, quiet movement and SilverSolutions in the powder coating to provide antimicrobial surface protection.

Basic American by Graham-Field is a leading manufacturer of beds serving the senior living industry for 40+ years. We proudly manufacture our beds in the USA, in Fond du Lac, Wisconsin, which enables faster response to customers' needs.

Innovation with industry-leading technology is designed into the Basic American Zenith, Matrix, and Liberty models to provide feature-rich beds and value for every budget.

6200
MATRIX BED
*shown in chair position with
optional two-pedal locking system*



MATRIX SERIES

FEATURES AND ADVANTAGES:

Made in the USA

- Minimal lead times
- Quality craftsmanship and manufacturing

SilverSolutions

- Located in the powder-coating which provides antimicrobial surface protection
- Premium Products

Sealed Ball Bearings

- Allows seamless, continuous, and noiseless movement
- Maintenance Free

Expandable Deck Available

- Pin-on wide kits extend the deck width to 39" or 42"

ANSI/AAMI ES60601 Safety Certification

- Tip stability testing
- Load capacity testing (twice the operating capacity)
- Roll clearance testing
- Splash proof electronics testing

Certified as a transport device for the CMS Emergency Preparedness Plan

- Easy maneuvering/adjustability through doorways and hallways during emergency or evacuation

Optional 'Plug and Play' digital electronics enable future feature additions

- Staff Control for Trendelenburg and Reverse Trendelenburg
- Chair Position
- Underbed Light
- USB Charger

Optional Intuitive two-pedal locking system

- Centrally located
- Easily visible and accessible at ALL heights and widths

Industry-leading, energy-saving electronics

- Consume up to 66% less power than transformer electronics

Grounded bed with right-angle hospital grade plug

- Provides added safety for both staff and residents

US Patented design of high/low mechanism

- Housed with the mainframe tubes; no holes are located on the top side that will harbor bacteria

Grid sleep surface with cross braces

- Sturdy, no-sag sleeping surface



6200 MATRIX BED
shown with pivot assist bar (front), counter-rotating assist (back) device available to add anytime at additional cost

USB power supply (right) underbed light (below)



EMBEDDED STAFF CONTROL
SHOWN ON 6200 MATRIX BED
available to add anytime at additional cost



6200 MATRIX BED
shown with pivot assist bar (front), counter-rotating assist (back) device

Specifications:

Adjustable Surface Length: <small>maximum</small>	84 inch
Adjustable Surface Length: <small>minimum</small>	80 inch
Fixed Sleep Deck Width:	36" or 42 inch
Height Travel Range: <small>maximum</small>	27 inch
Height Travel Range: <small>minimum</small>	8.95 inch
Safe Work Load: <small>includes resident, support surface, head/footboards, and accessories</small>	450 lb
Maximum Patient Weight: <small>EVENLY DISTRIBUTED</small>	400 lb
Actuators:	4 DC

Warranty:

Main Frame & Welds	15 years
Control Box & Actuator Motors	4 years
Hand Control Pendant & Cabling	3 years
Staff Control & Cabling	3 years
Head & Footboards	1 year
All Other Durable Components <small>not otherwise listed</small>	1 year

± 1 inch tolerance



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved. Basic American Medical Products is a trademark of GF Health Products, Inc.