

EXTENDED CARE BEDS



LIBERTY SERIES



The Liberty Series Bed provides you with the key features you want at a very affordable price point.

Basic American by Graham-Field is a leading manufacturer of beds serving the senior living industry for 40+ years. We proudly manufacture our beds in the USA, in Fond du Lac, Wisconsin, which enables faster response to customers' needs.

Innovation with industry-leading technology is designed into the Basic American Zenith, Matrix, and Liberty models to provide feature-rich beds and value for every budget.



LBT7174
LIBERTY BED

shown with optional head and foot end counter-rotating assists (sold separately)

LIBERTY SERIES

FEATURES AND ADVANTAGES:

Made in the USA

- Minimal lead times
- Quality craftsmanship and manufacturing

SilverSolutions

- Located in the powder-coating which provides antimicrobial surface protection
- Premium Products

Bed system UL rating

- Tip stability testing
- Roll clearance testing
- Splash proof electronics testing

Oversized Bearing Surfaces

- Reduces wear

Certified as a transport device for the CMS Emergency Preparedness Plan

- Easy maneuvering/adjustability through doorways and hallways during emergency or evacuation

Grounded bed with right-angle hospital grade plug


- Provides added safety for both staff and residents

Grid sleep surface with cross braces

- Sturdy, no-sag sleeping surface

Four Locking casters (standard)

- Optional pedal lock



Specifications:

Adjustable Surface Length:	80 inch
maximum	
Adjustable Surface Length:	76 inch
minimum	
Fixed Sleep Deck Width:	36 ¹ inch
Height Travel Range:	24.4 inch
maximum	
Height Travel Range:	9.95 inch
minimum	
Safe Work Load:	450 lb
includes resident, support surface, head/footboards, and accessories	
Maximum Patient Weight:	400 lb
EVENLY DISTRIBUTED	
Actuators:	3 DC

Warranty:

Main Frame & Welds	10 years
Control Box & Actuator Motors	2 years
Hand Control Pendant & Cabling	1 year
Staff Control & Cabling	1 year
Head & Footboards	1 year
All Other Durable Components	1 year
not otherwise listed	

± 1 inch tolerance



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website
www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Basic American Medical Products is a trademark of GF Health Products, Inc.

graham
field

GF2100057RevC25