



EXTENDED CARE BEDS

ZENITH SERIES

7200 & 9200 + SLIDE-W-I-D-E®



The Zenith™ Series by Basic American is the Premier Extended Care Bed. Engineered for durability and function, The Zenith features sealed ball bearings in the joints for smooth, quiet movement and SilverSolutions in the powder coating to provide antimicrobial surface protection.

Basic American by Graham-Field is a leading manufacturer of beds serving the senior living industry for 40+ years. We proudly manufacture our beds in the USA, in Fond du Lac, Wisconsin, which enables faster response to customers' needs.

Innovation with industry-leading technology is designed into the Basic American Zenith, Matrix, and Liberty models to provide feature-rich beds and value for every budget.

7200
ZENITH APS BED
*shown with pivot assist bar (right),
counter-rotating assist (left) device
available to add anytime at additional cost*

ZENITH SERIES

FEATURES AND ADVANTAGES:

Made in the USA

- Minimal lead times
- Quality craftsmanship and manufacturing

SilverSolutions

- Located in the powder-coating which provides antimicrobial surface protection
- Premium Products

Sealed Ball Bearings

- Allows seamless, continuous, and noiseless movement
- Maintenance Free

SLIDE-W-I-D-E® Integrated Expandable Deck

- Slides located on the frame effortlessly slide in and out allowing for quick and easy transition
- Expandable/retractable deck adjusts to 36", 39", and 42" with two quick-release ball-detent pins on each side

ANSI/AAMI ES60601 Safety Certification

- Tip stability testing
- Load capacity testing (twice the operating capacity)
- Roll clearance testing
- Splash proof electronics testing

Optional 'Plug and Play' digital electronics enable future feature additions

- Staff Control for Trendelenburg and Reverse Trendelenburg
- Chair Position
- Underbed Light
- USB Charger

Certified as a transport device for the CMS Emergency Preparedness Plan

- Easy maneuvering/adjustability through doorways and hallways during emergency or evacuation

Intuitive two-pedal locking system

- Centrally located
- Easily visible and accessible at ALL heights and widths

Industry-leading, energy-saving electronics

- Reduces Energy Consumption

Grounded bed with right-angle hospital grade plug

- Provides added safety for both staff and residents

US Patented design of high/low mechanism

- Housed with the mainframe tubes; no holes are located on the top side that will harbor bacteria

Grid sleep surface with cross braces

- Sturdy, no-sag sleeping surface



9200
ZENITH BED
shown with pivot assist bar (right),
counter-rotating assist (left) device
available to add anytime at additional cost



7200
ZENITH™ APS BED WITH
Slide-W-I-D-E
EXPANDABLE/RETRACTABLE SLEEP DECK

USB power supply (right)
underbed light (below)



Specifications:

	7200	9200
Adjustable Surface Length: maximum	84 inch	84 inch
Adjustable Surface Length: minimum	80 inch	80 inch
Fixed Sleep Deck Width:	36" inch	36" inch
Slide-W-I-D-E® option Sleep Deck Widths:	36" or 39" or 42 inch	36" or 39" or 42 inch
Height Travel Range: maximum	30 inch	30 inch
Height Travel Range: minimum	7.95 inch	7 inch
Safe Work Load: includes resident, support surface, head/footboards, and accessories	550 lb	600 lb
Maximum Patient Weight: EVENLY DISTRIBUTED	500 lb	550 lb
Actuators:	4 dc	4 dc

Warranty:

Main Frame & Welds	15 years	15 years
Control Box & Actuator Motors	5 years	5 years
Hand Control Pendant & Cabling	3 years	3 years
Staff Control & Cabling	3 years	3 years
Head & Footboards	1 year	1 year
All Other Durable Components not otherwise listed	1 year	1 year

± 1 inch tolerance



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

© 2021, GF Health Products, Inc. All Rights Reserved.
Basic American Medical Products and Slide-W-I-D-E are registered trademarks of GF Health Products, Inc.

graham
field

GF2100055Rev025