

RESPIRATORY

JB0112-070 NEBULIZER COMPRESSOR

NEB-U-TYKE® NEB-A-DOODLE™

This interactive, child-friendly design provides an ideal solution for effective aerosol treatment for children. Designed with an integrated drawing board, the Neb-a-Doodle will keep children relaxed and entertained during aerosol treatments and can be used multiple times. The Neb-u-Tyke® Neb-a Doodle Train is the latest addition to the Graham-Field line of quality nebulizer compressors that offer innovation, performance, and value for both patients and providers.

Features:

- Reliable, piston-driven compressor
- Interactive, drawing board for child to draw, doodle, and write on during aerosol treatment
- Convenient and compact design for easy transport and storage
- Comes complete with Nebulizer Accessory Kit
- Designed to meet the requirements of HCPCS Code: E0570
- Five year limited warranty



Package Contains:

- 1 Neb-u-Tyke® Neb-a-Doodle Nebulizer Compressor
- 1 User Manual
- 1 "T" Adaptor
- 1 Nebulizer Accessory Kit *complete with*
 - 1 Reusable nebulizer (*single patient use only*)
 - 1 Pediatric mask
 - 1 Adult Mask
 - 1 Mouthpiece
 - 1 PVC tubing
 - 1 Reservoir Tube
 - 3 Replacement filters



Options:

★ Item No. UOM Description

Item No.	UOM	Description
JB0112-070	1 ea	Neb-u-Tyke® Neb-a-Doodle Nebulizer
JB77400	1 ea	Neb-u-Tyke® Neb-a-Doodle & Case
JB76240	1 ea	Neb-a-Doodle Case ONLY

†Specifications:

Overall Width		Overall Length		Overall Height		Operating Range				Storage Range				Power	Product Weight	Sound Level	Max Pressure	Max Flow	Particle Size	Capacity						
low	high	low	high	low	high	low	high	low	high	low	high	low	high								compressor	nebulizer				
inch	cm	inch	cm	inch	cm	°F	°C	°F	°C	RH	RH	°F	°C	°F	°C	RH	RH	VAC	Hz	lb	kg	psi	lpm	µm	ml	
4.2	10.66	10	25.4	6.5	16.5	41	5	104	40	30%	95%	40	4.4	158	70	10%	95%	115	60	4	1.81	<57	32.5	≈12	≤5	10

★FDA (USA) restricts the sale of this device to or by the order of a physician

†Specifications are subject to change without notice.



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Neb-u-Tyke is a registered trademark of GF Health Products, Inc.