



READY SHIP PROGRAM

Basic American by Graham-Field offers several essential products within our Ready Ship Program. Ready Ship means the product you're looking for is readily available to make its way to you quickly and efficiently.

To help ensure that our Ready Ship initiative is seamless; certain ready ship items, such as dining chairs and casegoods are predesigned. This takes the guesswork out of the decision-making process and helps us deliver the needed goods in a timely manner.

Ready ship program includes: Beds and accessories, OBT's, Three Drawer Bedside, 4 Drawer Dressers, 24" Wardrobes, and Dining Chairs

Contact your sales rep for more details!



Crescent Cap headboard/
Crescent footboards shown
Custom orders only

READY SHIP BEDS & ACCESSORIES

We proudly manufacture all our beds right here in the USA! Innovation with industry-leading technology is designed into all Basic American Zenith, Matrix, and Liberty models to provide feature-rich beds.

Ready ship beds and accessories ship within 3-5 business days after order is released.

5 units max per order

Matrix 6200 SERIES

6200 WITH APS

ITEM N°: MZ6AF35QS

Crescent Cap headboard/Crescent footboards shown - custom orders only



Zenith 7200 SERIES

7200 ZENITH BED

ITEM N°: ZZ7F35QS

Crescent Cap headboard/Crescent footboards shown - custom orders only



Zenith 7200 SERIES

7200 SLIDE-W-I-D-E® WITH APS

EXPANDABLE/RETRACTABLE SLEEP DECK

ITEM N°: ZZ7ASWQS

Raised Panel head/footboards shown - custom orders only



LIBERTY BED

ITEM N°: LBTL7174QS

Center Crown head/foot boards shown



Board Options:



Rectangle Laminate Panel with Matching Edge




Center Crown Laminate Panel with Matching Edge


Finish Options:




Solar Oak



Cocobala



Wild Cherry



Windsor Mahogany

Bed Specifications:

	Zenith 7200	Matrix 6200	Liberty
Adjustable Patient Surface Length: maximum	84 inch	84 inch	80 inch
Adjustable Patient Surface Length: minimum	80 inch	80 inch	76 inch
Welded Sleep Deck Width:	35 inch	35 ^{or} 42 ^{inch}	35 inch
Height Travel Range: maximum	30 inch	27 inch	24.4 inch
Height Travel Range: minimum	7.95 inch	8.95 inch	9.95 inch
Safe Work Load: includes resident, support surface, head/footboards, and accessories	550 lb	450 lb	450 lb
Maximum Patient Weight: EVENLY DISTRIBUTED	500 lb	400 lb	400 lb
Actuators:	4 DC	4 DC	3 DC

Options & Features:

	Zenith 7200	Matrix 6200	Liberty
Chair Position (APSI)	■	■	—
Advanced Positioning Pendant: Reverse Trendelenburg/ Trendelenburg Positioning: via APS staff control 1	■	■	—
Slide-W-I-D-™ Option Available: between 35 - 42 inches upgrade controls anytime	■	—	—
Hand and/or Staff Controls: interchangeable, mix/match, standard/APSI	■	■	—
Head and Footboard Mounting:	Low Profile	Low Profile	Low Profile
Maximum Clearance: overbed tables and lifts	■	■	5.1 inch ⁴
FDA Guideline Compliance: for Bed Entrapment Risk Reductions	■	■	—
Silver Solutions Powder Coating: process provides antimicrobial surface protection ⁶	■	■	—
Made in the USA by Skilled Craftsmen: built to withstand trouble-free use for many years	■	■	—

Warranty:

	Zenith 7200	Matrix 6200	Liberty
Main Frame & Welds	15 years	15 years	10 years
Control Box & Actuator Motors	5 years	4 years	2 years
Hand Control Pendant & Cabling	3 years	3 years	1 year
Staff Control & Cabling	3 years	3 years	1 year
Head & Footboards	1 year	1 year	1 year

■ Standard ■ Available ■ Optional — Not Applicable/Available

Head & Footboards

Staff Control	7200 & 6200		inch	head & footboard finish			
	model	shape		Rectangle Center Crown Size	Solar Oak	Cocobala	Wild Cherry
Without	■	■	35	AB53600	AB20600	—	—
	■	■	42	MW53600	MW20600	—	—
Standard	■	■	35	ML53650	ML20650	—	—
	■	■	42	MWL20650	—	—	—
APS	■	■	35	MAL53650	MAL20650	—	—
	■	■	42	—	—	AB04600CCR	AB29600CCR
Without	■	■	35	—	—	MW04600CCR	MW29600CCR
	■	■	42	—	—	ML4650CCR	ML29650CCR
Standard	■	■	35	—	—	MAL4650CCR	MAL29650CCR
	■	■	42	—	—	MWAL04650CCR	MWAL29650CCR
Without	■	■	35	AB53600	AB20600	—	—
	■	■	35	—	—	AB04600CCR	AB29600CCR

■ Width Option ■ Board Style ■ Fits Bed — Not Available
ready ship head and footboards available to ship upon order release

Assist Devices:

	Zenith & Matrix	Liberty
	item no.	item no.
Counter-Rotating: ²	ZA91300	ZA82800
Pivoting: ³	ZA94700	ZA90350
Fixed: ³	ZA94900	ZA80300

Fixed Kit:


Zenith & Matrix
item no.
42" Wide Kit
ZA95600

(Converts 35" Fixed Width to 42" Fixed Width)

Counter-Rotating Assist Device



Fixed Assist Devices



Pivot Assist Bars



Crescent Cap headboard shown - custom orders only

¹ anytime add-on or with purchase ² counter rotating assist device for head end - uow /pair ³ both pivot & fixed assist devices are tool-less ⁴ fixed base ⁵ Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment™ details: www.fda.gov ⁶ preserves the cured coating during the service life of the bed & inhibits the growth of organisms that can detract from the appearance and coating functionality, inhibits the growth of bacterial odors and resists stains caused by mildew and/or mold



READY SHIP OBT


We offer ready ship 3DL overbed tops in four of our most popular finishes, offering both a raised edge and flat edge option. Along with the top, you'll receive a standard 40lb base or 100lb* base depending on what is ordered.

Ready ship OBT's ship within 3-5 business days after order approval.
5 units max per order



*100 lb base only compatible with oversized 18" x 35"

AVAILABLE FINISHES




Cocobala



Solar Oak



Wild Cherry



Windsor Mahogany

	Item No.	A793053	A793020	A793029	A793004	A860053	A860029	A860004	A797053	A797020	A797029	A797004	A859053	A859029	A859004	Oversized Tops			
		A919053XL	A919020XL	A919029XL	A919004XL														
Low Height ¹	base height	■	■	■	■	■	■	■	—	—	—	—	—	—	—	—	—	—	—
Standard Height ¹	base height	—	—	—	—	—	—	—	■	■	■	■	■	■	■	—	—	—	—
Extended Height ¹	base height	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—
H Base	base style	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—
Modified H Base	base style	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
material	3DL	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Flat Edge	top profile	■	■	■	■	—	—	—	■	■	■	■	—	—	—	—	—	—	—
Raised Edge	top profile	—	—	—	—	■	■	■	—	—	—	—	■	■	■	—	—	—	—
Solar Oak	top finish	■	—	—	—	■	—	—	■	—	—	—	■	—	—	—	—	—	—
Cocobala	top finish	—	■	—	—	—	—	—	—	■	—	—	—	—	—	—	—	—	—
Windsor Mahogany	top finish	—	—	■	—	—	■	—	—	—	■	—	—	■	—	—	—	—	—
Wild Cherry	top finish	—	—	—	■	—	—	■	—	—	—	■	—	—	■	—	—	—	■
Weight Capacity EVENLY DISTRIBUTED	lb	40	40	40	40	40	40	40	40	40	40	40	40	40	40	100	100	100	100
UOM	ea	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1

¹ height adjusts infinitely

■ Standard Feature ■ Available Style ■ Available Finish — Not Available/Does Not Apply



READY SHIP CASEGOODS

We are happy to offer the Trinity Park casegood collection within our ready ship program. This group offers simplistic styling while also providing a practical solution to resident room furnishing.


3DL Option: Trinity Park - 15 units max per order in any combination. Units ship within 15 business days after order is released. This collection is backed by an 8-year warranty.

Available Units:


Trinity Park:

3 Drawer Bedside	A31-37-RS
4 Drawer Dresser	A37-37-RS
24" Two Door/Two Drawer Wardrobe	B02-371-RS
Hardware	Bolt Strap Pull Brushed Silver


AVAILABLE FINISHES



Talas Cherry



Shaker Cherry



Cocobala



*BOLT STRAP HARDWARE ONLY
AVAILABLE ON TRINITY PARK



READY SHIP CHAIRS

We are thrilled to provide dining chair options that are predesigned and ready to ship! Our lightweight, aluminum frames are our best sellers and for good reason. These frames have the “look” of a wood frame chair while providing unmatched durability and are backed by a 12-year frame warranty.

Baxter and Ryno - 20 units max per order. Chairs ship within 3 business days after the order is released. Fabric and finish on chairs are sold/shipped as shown.

Available Styles:

	Item:	Finish:
Baxter	249-F10-QS4	Dark Walnut
Ryno	471-98-QS1-1	Slate




BAXTER




RYNO



FINISHES



Dark Walnut
*Baxter



Slate
*Ryno

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

1.770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2023, GF Health Products, Inc. All Rights Reserved.
Basic American Medical Products is a trademark of GF Health Products, Inc.

grahamfield

GF2300025RevE26