

# RESPIRATORY

## JB0112-071 NEBULIZER COMPRESSOR NEB-U-TYKE® BELLA BUTTERFLY

This interactive, child-friendly design provides an ideal solution for effective aerosol treatment for children. Designed with an attractive butterfly shape and colorful stickers to personalize, Bella Butterfly will keep child relaxed and entertained during aerosol treatment. The Neb-u-Tyke® Bella Butterfly is the latest addition to the Graham-Field line of quality nebulizer compressors that offer innovation, performance, and value for both patients and providers.

### Features:

- Reliable, piston-driven compressor
- Child-friendly design complete with stickers to personalize during aerosol treatment
- Convenient and compact design for easy transport and storage
- Piston Driven Method of Nebulization
- Comes complete with Nebulizer Accessory Kit
- Designed to meet the requirements of HCPCS Code: E0570
- Five year limited warranty



### Package Contains:

- 1 Neb-u-Tyke® Bella Butterfly Nebulizer Compressor
- 1 User Manual
- 1 Nebulizer Accessory Kit *complete with*
  - 1 Reusable nebulizer (*single patient use only*)
  - 1 Pediatric mask
  - 1 Mouthpiece
  - 1 PVC tubing
  - 5 Replacement filters



### Options:

★ Item No.	UOM	Description
JB0112-071	1 ea	Neb-u-Tyke® Nebulizer
JB0112-071F	5 /pk	Replacement Filters

### †Specifications:

Overall Width		Overall Length		Overall Height		Operating Range				Storage Range				Power		Product Weight		Sound Level	Max Pressure	Max Flow	Particle Size Capacity			
inch	cm	inch	cm	inch	cm	Temperature		Humidity		Temperature		Humidity		VAC	Hz	lb	kg	dBA	compressor	nebulizer				
						low	high	low	high	low	high	low	high						psi	lpm	µm	ml		
6.34	16.1	9.10	23.1	3.70	9.3	50 °F	104 °F	10% RH	95% RH	4 °F	-15 °F	140 °F	60 °F	10% RH	95% RH	120	60	2.88	1.31	<60	42	≈8	≤5	5

★FDA (USA) restricts the sale of this device to or by the order of a physician

†Specifications are subject to change without notice.

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. [www.grahamfield.com](http://www.grahamfield.com)



© 2013, GF Health Products, Inc. All Rights Reserved. Neb-u-Tyke is a registered trademark of GF Health Products, Inc.