


# REPLACEMENT TOPS



## REPLACEMENT TOPS

Make your exam tables look new at a fraction of the cost. All components within the replacement top are all new including wood, foam, vinyl and hardware.

- Industry standard colors *plus*
- Replacement Tops fit leading brand tables




**GF100RT-xx**  
shown in violet | VI  
**100**




**GF104RT-xx**  
shown in teal | TE  
**104**



**GF204RT-xx**  
shown in british tan | BT  
**204**



**GF304RT-xx**  
shown in ocean gray | OG  
**304**



**GF404RT-xx**  
shown in doe | DO  
**404**



**GF405RT-xx**  
shown in shadow | SH  
**405**



**GF420RT-xx**  
shown in artichoke | AR  
**420**



**GF460RT-xx**  
shown in british tan | BT  
**460**

### Specifications:

	100   104		204 <sup>1</sup>		304		404   405		420   460	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
Overall Width:	27.50	69.85	27.50	69.85	27.50	69.85	27.50	69.85	27.50	69.85
Overall Length:	55.00	139.70	55.00	139.70	55.00	139.70	60.00	152.40	60.00	152.40
Footrest Width:	20.00	50.80	19.00	48.26	20.00	50.80	20.50	52.07	19.00	48.26
Footrest Length:	18.00	45.72	18.00	45.72	18.00	45.72	17.50	44.45	18.00	45.72
Shipping Weight:	lb	kg	lb	kg	lb	kg	lb	kg	lb	kg
	55	24.95	55	24.95	55	24.95	35	15.88	55	24.95
<b>Features:</b>										
Hook & Loop Attachment:	—		—		—		included		—	
Paper Roll Holder:	included		—		—		—		included	
Mounting Hardware:	included		—		—		—		included	
Lift Handles:	—		not included <sup>1</sup>		—		—		—	
Top Replacement Knobs:	—		not included <sup>1</sup>		—		—		—	
Individual Product Links:	<a href="#">🔗</a> <a href="#">🔗</a>		<a href="#">🔗</a>		<a href="#">🔗</a>		<a href="#">🔗</a> <a href="#">🔗</a>		<a href="#">🔗</a> <a href="#">🔗</a>	

<sup>1</sup> 2009 year model and later do not include lift handles and top attachment knobs  
If the OEM model has a seamless top, request paper roll holder when ordering  
— Does Not Apply    🔗 Active Hyperlink

### Standard Color Options (-xx)

Available in 24 color choices to complement most decors.

AL ALABASTER	AQ AQUA GREEN	AR ARTICHOKE	BL BLACK	BU BLUE	BT BRITISH TAN
DO DOE	DB DUSTY BLUE	GR GRAY	MO MOSS	NA NAVY	NO NORTHWOODS
OG OCEAN GRAY	PM PALM	PL PEARL	PU PLUM	ST SADDLE TAN	SH SHADOW
WH SUPER WHITE	TL TEA LEAF	TE TEAL	VI VIOLET	WG WARM GRAY	WW WEDGEWOOD

Swatches are digital reproductions. perceived texture and coloring may vary on finished product. Monitor and printer settings will effect appearance of color. Swatch Card is available from your representative.

GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.hausted.com](http://www.hausted.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2018, GF Health Products, Inc. All Rights Reserved.  
Hausted is a registered trademark of GF Health Products, Inc.