

ULTRA-LOW SLEEP SYSTEM

AMERILUX™

Introducing AmeriLux™ — Where Safety Meets Style

We're excited to unveil AmeriLux™, the latest innovation in our extended care bed line—designed to elevate resident safety while blending seamlessly into residential-style settings.

With its ultra-low profile and sleek black frame, AmeriLux™ redefines the look of traditional beds—offering a more modern, design-friendly aesthetic without compromising on safety.

Built with advanced Ultra-Low technology, AmeriLux™ significantly reduces the risk of fall-related injuries by lowering close to the ground—providing peace of mind for both caregivers and families.

For residents at higher risk, an optional fall mat transforms potential falls into gentle rolls, offering an added layer of protection without feeling clinical.

With AmeriLux™, you don't have to choose between elegance and safety—you get both.





AMERILUX SERIES

FEATURES AND ADVANTAGES:

Made in the USA

- Minimal lead times
- Quality craftsmanship and manufacturing

Ultra-Low Technology

- Adjustable height from floor level at 4 inches to nursing height at 30 inches
- Bed electronics feature full floor clearance, preventing floor contact vs similar models offered by competitors
- Available as 36 inch width or 42 inch width bed
- Sleek, black frame
- Advanced Positioning System (APS) includes One-Button Comfort/Cardiac Chair Position (Embedded staff control in footboard allows Trendelenburg/Reverse Trendelenburg)

Customizable Design Options

- Creative flexibility to align seamlessly with your design aesthetic
- Personalization by adding wooden side boards for a cozy, home-inspired feel
- Offered in various board styles & multiplies finishes, with the possibility of custom styles or finishes upon request

SilverSolutions™

- Located in the powder-coating which provides antimicrobial surface protection
- Deters the growth of odor-causing bacteria

Sealed Ball Bearings

- Allows seamless, continuous, and noiseless movement
- Maintenance Free

Optional 'Plug and Play' digital electronics enable future feature additions

- Staff Control for Trendelenburg and Reverse Trendelenburg
- Under-Bed Light
- USB Charger

Caster locking system

- Four locking casters (standard)
- Easily visible and accessible at ALL heights and widths

Industry-leading, energy-saving electronics

- Reduces Energy Consumption

Grounded bed with right-angle hospital grade plug


- Provides added safety for both staff and residents

US Patented design of high/low mechanism

- Housed with the mainframe tubes; no holes are located on the top side that will harbor bacteria

Grid sleep surface with cross braces

- Sturdy, no-sag sleeping surface



Fall Mat

USB Power Supply

Under-Bed Light

Fixed Assist



ALX36
AMERILUX BED
shown with counter-rotating assist device

Specifications:

Adjustable Surface Length: <small>maximum</small>	84 inch
Adjustable Surface Length: <small>minimum</small>	80 inch
Fixed Sleep Deck Width:	36" or 42 inch
Height Travel Range: <small>maximum</small>	30 inch
Height Travel Range: <small>minimum</small>	4" inch
Safe Work Load: <small>includes resident, support surface, head/footboards, and accessories.</small>	550 lb
Maximum Patient Weight: <small>EVENLY DISTRIBUTED</small>	500 lb
Actuators:	4 DC

Warranty:

Main Frame & Welds	12 years
Control Box & Actuator Motors	5 years
Hand Control Pendant & Cabling	3 years
Staff Control & Cabling	3 years
Head & Footboards	1 year
All Other Durable Components <small>not otherwise listed</small>	1 year

± 1 inch tolerance



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com



© 2025, GF Health Products, Inc. All Rights Reserved.
Basic American Medical Products is a trademark of GF Health Products, Inc.