

EXAM ROOM SUPPLIES

# MEDIUM DUTY CART

This versatile medium duty cart with side rails easily moves records, supplies and equipment throughout your facility. It's made of 16GA stainless steel, has three solid shelves with lipped edges on all four sides, travels on four 4" swivel casters and has a weight capacity of 500 lb evenly distributed.

### Features

- Stainless steel construction for easy cleaning
- Three Shelves
- One Push handle
- 4" (10.2 cm) swivel casters



	overall Width		overall Height		overall Depth		Caster diameter		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	
GF-2427535000	27.5	69.9	40	101.6	18.125	46	4	10.2	500	226.8	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2025, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.