

EXAM ROOM SUPPLIES

BENJAMIN MAYO STANDS

The Benjamin Mayo stand design is number one in ease and durability. Its state-of-the-art double post design minimizes lateral movement and increases overall unit stability. Via a foot-operated mechanism, you can adjust the height from 39 ½" to 62" high off the ground. A guard over the foot prevents accidental release.

Features

- Double post stainless steel Mayo stand
- Double post design minimizes lateral movement and increases overall stability
- Two-pronged base with two 2" (5.1 cm) roller blade back casters
- Three tray sizes:
 - Standard: 12 5/8" x 19 1/8"
 - Medium: 16" x 21"
 - Large: 20" x 25"
- Foot-operated lift mechanism with a guard to prevent accidental release
- Height Adjustment: 39 ½" - 62"



	Tray				Adjustable Height				Weight Capacity		UOM
	width		depth		minimum		maximum		EVENLY DISTRIBUTED		
	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	
GF-0668867000	12.625	32.1	19.125	48.6	39.5	100.3	59.5	151.1			1 ea
GF-0668869000	16	40.6	21	53.3	39.5	100.3	59.5	151.1			1 ea
GF-0668871000	20	50.8	25	63.5	39.5	100.3	59.5	151.1			1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2025, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.