

EXAM ROOM SUPPLIES

LIGHT DUTY CART

This versatile light duty cart easily moves records, supplies and equipment throughout your facility. It's made of 20GA stainless steel, has three solid shelves, one push handle with protective rubber bumpers on each corner, travels on four 4" swivel casters and has a weight capacity of 300 lbs. evenly distributed. The top and middle shelves have three edges with 1¼" sides and a front open edge to slide equipment in and out with ease. The bottom shelf has four open edges for greater accessibility to equipment stored in a lower position.

Features

- Stainless steel construction for easy cleaning
- Three Shelves
- One Push handle
- 4" (10.2 cm) swivel casters



	overall Width		overall Height		overall Depth		Caster diameter		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	
GF-2427534002	28.25	71.8	35.75	90.8	16	40.6	4	10.2	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2025, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.