



## HEALTHCARE SEATING

### 6200 SERIES

# PHLEBOTOMY CHAIR

The Lumex 6200 Series bariatric blood draw/phlebotomy chair is designed with both the patient and phlebotomist in mind. The sensible and affordable design offers a small footprint and modern style while keeping patient comfort and well-being the highest priority.

#### Features:

- Fully upholstered padded seat, back, and arm rests
- Individual height adjustable arm rests
- Wallsaver wireframe included
- Powder coat frame finish - Foghorn
- Available in 24 industry standard colors to complement most decors

[6200 Product Info](#)

[Fabric Choices](#)



#### 6200 Options:

6200-xx	<input checked="" type="checkbox"/>	Single Flip Arm
6200DA-xx	<input checked="" type="checkbox"/>	Dual Flip Arm

★ options available at additional cost



#### Specifications:

Overall Width		Overall Depth		Overall Height		Width between arms		Arm Width		Arm Height <sup>8</sup>				Seat Width		Seat Depth		Seat Height		Back Width		Back Height from seat		Product Weight <sup>8</sup>		Shipping Weight <sup>8</sup>		Weight <sup>1</sup> Capacity	
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	lb	kg†	lb	kg†	lb	kg		
35.00	88.90	28.00	71.12	37.50	95.25	25.00	63.50	5.00	12.70	26.50	67.31	32.00	81.28	25.50	64.77	18.75	47.63	22.00	55.88	26.00	66.04	16.00	40.64	52	23.59	55	24.95	800	362.87

★ options available at additional cost †-xx in the product number is the placeholder for fabric choice § images show this accessory in use - sold separately at additional cost  
 † standard rounding method used ‡ maximum weight capacity - evenly distributed - refer to user manual § approximate measurement - options may alter spec



GF Health Products, Inc. is an  
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
 The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
 Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF1800043 Rev1/26