

EXTENDED CARE BED ACCESSORIES

HIGH FALL PAD

The **Basic American High Fall Pad** is designed for use with the **Basic American AmeriLux bed** in its lowest height position. Acting as an essential accessory, the pad provides an extended surface at about the same height as the bed and a 6" standard mattress.

This design ensures that in the event of a fall, patients can safely roll onto the pad, reducing the risk of injury while offering peace of mind. By enhancing safety, the pad may also contribute to fewer sleep disruptions, helping individuals feel more secure throughout the night.

Key Features

- Extended surface at mattress height for seamless protection
- Roll-off design can reduce injury risk
- Enhances patient comfort and security
- Compatible with the AmeriLux bed
- Built-in handles for easy maneuvering and storage transport
- Buckle closure ensures the pad stays compactly folded for storage and transport

Warranty

- **Foam:** Five (5) years
- **Cover:** Two (2) years



High Fall Pad is not intended to be used as a sleeping surface.

Measurements & Specifications

	overall Width		overall Height		overall Length		Product Weight		UOM
	inch	cm	inch	cm	inch	cm	lb	kg	
FP3076	30	75	9	22.8	76	193	29	13.2	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2026, GF Health Products, Inc. All Rights Reserved.
Basic American Medical Products is a trademark of GF Health Products, Inc.