



WHEELCHAIRS

STEEL TRANSPORT CHAIR

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device. Its compact size and lighter weight allow a companion to lift and store the chair with ease.

Features:

- Available in 17" and 19" seat widths
- Steel frame in hammertone silvervein finish
- Fold-down back makes transport and storage easy
- Black padded nylon upholstery is flame retardant
- Carry straps on seat allow for ease of lifting and folding
- Full length padded arms provide user comfort
- Tool-free, adjustable, detachable swingaway footrests
- Black padded nylon 3" leg strap standard. Velcro adjustable loops slide over installed footrests
- Attendant wheel locks on both rear casters for better control and access by companion
- One year limited warranty

Product Info



Item	HCPCS ²	Overall Chair Width		Overall Chair Length		Overall Chair Height		Folded Width		Seat Width		Seat Depth		Seat Height from floor		Caster Diameter		Shipping Weight ⁷		Product Weight ⁸ w/o front rigging		Weight Capacity ¹ EVENLY DISTRIBUTED	
		inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	lb [†]	kg [†]	lb [†]	kg [†]	lb	kg [†]
EJ795-1	E1038	22.00	55.88	37.00	93.98	39.00	99.06	10.00	25.40	17.00	43.18	15.75	40.00	19.00	48.26	8.00	20.32	33	14.97	24	10.89	250	113.40
EJ796-1		24.00	60.96							19.00	48.26									29	13.15		

all dimensions ±0.25" variance (0.64cm) † standard rounding method used ¹ maximum weight capacity - evenly distributed - refer to user manual ³ designed to meet the requirements of HCPCS code ⁸ approximate measurement - size affects spec

Parts & Accessories: ★ sold separately

	HCPCS ³	UOM	Item No.
42" Positioning Belt (plastic buckle)	E0978	1/ea	EJ76-BELT



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

+1 770.368.4700

Information contained herein is subject to change without notice. The most current & complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.
Velcro is a registered trademark of Velcro IP Holdings LLC

graham field™

GF1000034Rev G26