

PATIENT LIFTS

BARIATRIC TRAPEZE

Helps patient safely change position in bed and transfer from bed to other equipment with minimal assistance.

Features:

- Rubber bumper guards protect walls
- Handbar location fully adjustable
- Comes with floor stand
- 14mm wrench required for assembly

 Product Info



Item	HCPCS ³	Overall Width		Overall Length		Overall Height		Product Weight ⁸		assisted Weight Capacity ¹	
		inch	cm†	inch	cm†	inch	cm†	lb	kg†	lb	kg†
2940B	E0912	28.88	73.35	40.50	102.87	76.88	195.28	72	32.66	450	204.12
2960B	E0912	29.25	74.30	41.00	104.14	80.63	204.80	90	40.82	600	272.16

¹ maximum weight capacity - evenly distributed - refer to user manual ³ designed to meet the requirements of HCPCS code x
⁸ approximate measurement - options affect spec [†] standard rounding method used



GF Health Products, Inc. is an
 intertek ISO 13485:2016 / MDSAP Certified Company

Images are digital reproductions. perceived texture and coloring may vary on finished product. Monitor and printer settings will effect appearance of color.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current & complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
 Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0700060RevC26