

PATIENT LIFTS

FLOOR STAND

Provides freestanding capability for Versa-Helper Trapeze 2800x

Features:

- Individually-leveled floor bars provide firm, secure base
- Folds flat for easy storage and handling
- Floor Stand for use with Versa-Helper Trapeze designed to meet the requirements of HCPCS code: E0940

Product Info

2800 Compatible Trapeze Info



Item	Finish	HCPCS ³	Overall Width		Overall Length		Overall Height		Width between legs		assisted Weight ¹ Capacity	
			inch	cm†	inch	cm†	inch	cm†	inch	cm†	lb	kg†
2840A	Chrome-Plated	E0940	29.20	74.17	40.00	101.60	39.00	99.06	27.00	68.58	250	113.34
2840GA	Gray	E0940	29.20	74.17	40.00	101.60	39.00	99.06	27.00	68.58	250	113.34

¹ maximum weight capacity - evenly distributed - refer to user manual

³ designed to meet the requirements of HCPCS code when used with Versa-Helper Trapeze 2800x

⁸ approximate measurement - options affect spec

[†] standard rounding method used



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

+1 770.368.4700

Information contained herein is subject to change without notice. The most current & complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0700057/RevC26