

MOBILITY



CANES

The Lumex® Canes are stylish, unique and are available in a wide assortment of handgrip styles and sizes to fit a broad range of users with ease. The canes are offered in both standard and bariatric weight capacities and designed for both right and left-hand use.

WOOD CANES

- Crafted in satin-smooth hardwood with durable rubber tip
- Available in Standard and Derby Style
- Item No. 5180A, 5182A, 5183A, 5184A

ALUMINUM CANES

- Contoured vinyl or foam handle provides a comfortable, secure handhold
- Attractive finish resists marring and scratching
- Easy push-button height adjustment
- Locking collet keeps tubing secure
- Item No. 6220A, 6221A, 6223A, 6224A

FOLDING CANES

- Folds to 11' length for convenient transport
- Derby-style handle is contoured for comfort
- Easy push-button height adjustment
- Comes with convenient carry pouch
- Item No. 5950A-4

TRI-STEP FOLDING CANE

6000 Series

- Three point contact base allows cane to stand upright
- Comfortable, contoured grip
- Lightweight, aluminum folding design



E0100*



E0100*



E0100*




E0100*

OFFSET CANES

- Offset handle style is designed with handle that places user's weight directly over the base for maximum balance and control
- Attractive finish resists marring and scratching
- Easy push-button height adjustment
- Locking collet keeps tubing secure
- Ortho-Ease grip is ideal for arthritis and those needing a wider grip for support (6327 & 6328)
- Item No. 5940A, 5941A, 5942A, 5943A, 5944A, 5950BI, 6327, 6328, 5950SB, 5950SBK, 5950SR



E0100*



E01100*

AMERICAN PRIDE OFFSET CANES

- Show American and military pride spirit
- Soft, foam grip provides support
- Item No. 5970CA & 5970FL

QUAD CANES

- Low profile, four-point base is designed for maximum stability and walking assistance
- Rounded-edge base helps prevent bruising and scraping of ankles
- Easy push-button height adjustment
- Locking collet keeps tubing secure
- Can be adjusted to use properly for left or right-handed user
- Item No. 6140A, 6140BK, 6141A, 6141BK, 6120A, 6121A, 6125A, 6130A (Bariatric)



E0105*

*HCPCS Code

Item No.	Adjustment Range	Finish	Handgrip Style	Weight Capacity ¹ EVENLY DISTRIBUTED	UOM
5180A, 5182A, 5183A	36 in / 91.4 cm	Walnut, Mahogany	Standard	250 lb / 113.3 kg	6 /cs
5184A	42 in / 106.6 cm	Walnut	Standard	250 lb / 113.3 kg	6 /cs
5950A-4	32.5 in - 36.5 in 82.5 cm - 92.7 cm	Black	Derby	250 lb / 113.3 kg	4 /cs
5940A, 5941A, 5942A, 5943A, 5944A	31 in - 39 in 78.7 cm - 99 cm	5940A/5943A - Silver; 5941A, 5942A - Black; 5944A-Multipack;	5940A: PVC; Other models: Nitrile Foam	250 lb / 113.3 kg	6 /cs
5950BI, 5950SB, 5950SBK, 5950SR, 5970CA, 5970FL	30 in - 39 in 76.2 cm - 99 cm	5950BI - Blue Ice, SB (Blue), BK (Black), SR (Red), SS (Silver), Flag, Camo	Nitrile Foam	250 lb / 113.3 kg	6 /cs
6327, 6328	34 in - 42 in 86.3 cm - 106.6 cm	Bronze	Nitrile Foam	250 lb / 113.3 kg	6 /cs
6334A	29 in - 38 in 73.6 cm - 96.5 cm	Black	Nitrile Foam	600 lb / 272.1 kg	6334A-1 : 1 /ea 6334A-4 : 4 /cs
6140A, 6140BK, 6141A, 6141BK, 6120A, 6121A, 6125A	29 in - 38 in 73.6 cm - 96.5 cm	6140A, 6141A, 6120A, 6121A - Aluminum; 6140BK, 6141BK - Black 6125A - Bronze	Standard PVC (6120A, 6121A, 6124A, 6125A - Ortho-Ease)	300 lb / 136 kg	4 /cs
6130A	29 in - 38 in 73.6 cm - 96.5 cm	Silver Black Hammertone	Standard PVC	500 lb / 226.7 kg	2 /cs
6000B, BK, G, P	33 in - 37 in 83.8 cm - 93.9 cm	Blue, Black, Green, Pink	Derby	300 lb / 136 kg	4 /cs

CRUTCHES

The Lumex Crutches offer a solid one-piece crutch frame in lightweight anodized aluminum with push-button height adjustment and contoured handgrips. The handgrip position is readily adjusted with wingnuts and has the push-button height adjustment and patient height guide to allow the crutch to be fitted to the proper user height. The Lumex Forearm Crutches are constructed of high-strength aluminum tubing, vinyl coated arm cuffs that can be shaped to the user's arms, cushioned vinyl handgrips and extra large crutch tips that provide added traction and stability.

UNIVERSAL ALUMINUM CRUTCHES

- Reinforced center tube and footpiece provide added strength
- Cushioned underarm pad and handgrips for user comfort
- Easy to adjust handgrip with wingnut to secure positioning
- Item No.
3610LF-8, 3611LF-8, 3612LF-8, 3613C-8



E0114*

DELUXE ORTHO FOREARM CRUTCHES

- Offers all of the features of the Deluxe Forearm crutches with the added benefit of the therapist-preferred Ortho-Ease handgrip
- Handgrip angle can be adjusted to comfortably fit the user's hand and wrist angle
- Handgrip adjustment can properly position the handgrip to provide the correct elbow flexion which can reduce fatigue and minimize stress on the joints
- Available in Small, Medium and Large sizes
- Item No. w, 6346, 6347



E0110*

Product Webpage Inactive

DELUXE FOREARM CRUTCHES

- Designed with vinyl coated arm cuffs that can be shaped to user's arms for a custom fit
- Leg and forearm sections independently adjust with easy push-button adjustment in 1" increments
- Extra large crutch tips provide added traction and stability
- Contoured vinyl handgrips for user comfort
- Available in Adult, Tall and Youth sizes
- Item No. 6350A, 6350T, 6350Y



E0110*



*HCPCS Code

Item No.	Arm Height Adjustment	Armpad/ Cuff to Floor Height Adjustment	Recommended Patient Height	Weight Capacity ¹ EVENLY DISTRIBUTED	UOM
3610LF-8 (Adult)	15"-19" (Adult)	45"-53" (A)	5' 2"-5' 10" (A)	300 lb / 136 kg	8 pr / cs
3611LF-8 (Tall)	16"-20" (Tall)	53"-61.5" (T)	5' 10"-6' 6" (A)		
3612LF-8 (Youth)	9"-13" (Youth)	37"-46" (Y)	4' 6"-5' 2" (Y)		
3613C-8 (Combo)					
6350A (Adult)	9" - 12"	42"-51" (T)	5' 10"-6' 6" (T)	300 lb / 136 kg	1 pr
6350T (Tall)	from handle to top of cuff	37.5"-46.5" (A)	5' 0"-6' 2" (A)		
6350Y (Youth)		31"-40" (Y)	4' 2"-5' 2" (Y)		
6345 (Large)	9.2" - 12.2"	42"-51" (L)	5' 10"-6' 6" (L)	300 lb / 136 kg	1 pr
6346 (Medium)	from handle to top of cuff	37.5"-46.5" (M)	5' 0"-6' 2" (M)		
6347 (Small)		31"-40" (S)	4' 2"-5' 2" (S)		

¹ maximum weight capacity - evenly distributed - refer to user manual

Lumex Mobility Accessories



Walker Pouch



Walker Glide Skis



Walker Tray



Mobility Lights

The Lumex Folding Walkers are designed to offer consumers and clinicians a high quality walking aid that gives the user more independence in their everyday lives. The Lumex walkers are constructed using lightweight aluminum making them easy to lift and maneuver. The walkers feature side braces for stability and a secure handhold and are offered in both adult and junior sizes.

Dual-Release Folding Walker

- Raised H-frame design facilitates over-toilet positioning
- Rigid dual side brace design
- Dual-release frame locking mechanism
- PVC handgrips for comfort and long lasting wear
- Easy push-button height adjustment on leg attachments
- Available with or without 5" front fixed wheels
- Rear glide caps to facilitate use on most surfaces for wheeled options
- Item No. 616070A, 616270A-4, 616370A-4



E0135*

Single-Release Folding Walker

- Single release folding mechanism is designed to aid users with limited hand dexterity
- Raised H-frame design facilitates over-toilet positioning
- Rigid lower side brace design
- PVC handgrips for comfort and long lasting wear
- Easy push-button height adjustment on leg attachments
- Item No. 613070A



E0135*

Everyday Dual-Release Folding Walker

- Dual-release frame locking mechanism
- PVC handgrips for comfort and long lasting wear
- Easy push-button height adjustment on leg attachments
- Available with or without 5" front fixed wheels
- Rear glide tips to enhance maneuverability on most surfaces for wheeled options
- Item No. 716070A-2, 716070A-4, 716170A-4, 716270A-2, 716270A-4, 716370A-2, 716370A-4



E0135*

Rigid Walker

- Non-folding design provides extra rigidity for use where space is not an issue
- Angled adjustable legs provide maximum patient support in all directions and help resist tipping
- Angled side brace helps support patients as they stand up and sit down
- Item No. 1925A



E0135*

Bariatric Dual-Release X-Wide Walker

- Increased width and depth accommodates larger users
- Frame is designed to be robust and rigid to provide additional strength
- Available with and without double 5" fixed wheels
- Item No. 604070A, 604070W (with wheels)



E0135*

colorSelect

ColorSelect Dual-Release Folding Walker

- Same features as the Everyday Dual-Release Walker but available in a variety of color options to appeal to the user
- Item No.* 716270B-2, 716270BK-2, 716270G-2, 716270R-2, 716270P-2, 716270PK-2
- *1 ea also available



Platform Walker Attachment

- Designed for users with limited strength or dexterity who are unable to grip the walker
- Adjustable strap comfortably secures arm in the forearm rest
- Padded forearm rest can be pivoted and angled to fit user comfortably and to assist in lifting walker
- Universal clamp fits most 1" tubing walkers
- Tool-free installation and adjustment allows quick fitting for most patient sizes and types
- Item No. 6132B



UpRise Onyx Folding Walker

- Combination of folding walker, rising aid and toilet frame all in one
- Secondary handles provide stable assistance from a seated to a standing position
- Attractive, upscale Onyx finish
- Single release folding mechanism
- Durable aluminum tubing provides strength while remaining lightweight
- Adjustable height fits a broad range of user heights
- Item No. 700175C-2, 700175CR (Retail Packaged)



Walkane Hemi Walker

- Lighter than a walker, more stable than a quad cane
- Ideal for users who do not require a walker but need a wider base for support and for individuals with the use of only one hand or arm
- Contoured, foam hand grips are comfortable and provide secure hand hold
- Comfortable, dual level handgrips provide support from a seated to standing position
- Item No. 6015A



*HCPCS Code

Item No.	Width Inside Handgrips	Overall Width	Width Inside Rear Legs	Overall Depth	Minimum Height	Maximum Height	Weight Capacity ¹ EVENLY DISTRIBUTED
616070A	16.5 in / 41.9 cm	22.7 in / 41.9 cm	20.1 in / 51 cm	16.5 in / 41.9 cm	32.5 in / 82.5 cm	41.5 in / 105.4 cm	300 lb / 136 kg
616270A-4	16.5 in / 41.9 cm	26 in / 66 cm	20.1 in / 51 cm	16.5 in / 41.9 cm	34.5 in / 87.6 cm	41.5 in / 105.4 cm	300 lb / 136 kg
616370A-4	16.5 in / 41.9 cm	25.3 in / 64.2 cm	18.6 in / 47.2 cm	15.5 in / 39.3 cm	29.5 in / 74.9 cm	34.5 in / 87.6 cm	300 lb / 136 kg
604070A	20 in / 50.8 cm	26.1 in / 66.2 cm	23 in / 55.4 cm	18 in / 45.7 cm	31 in / 78.7 cm	38 in / 96.5 cm	600 lb / 272.1 kg
604070W	20 in / 50.8 cm	29.2 in / 74.1 cm	23 in / 55.4 cm	19 in / 48.2 cm	33 in / 83.8 cm	40 in / 101.6 cm	600 lb / 272.1 kg
716070A-1 716070A-2 716070A-4	17 in / 43.1 cm	22 in / 55.8 cm	20 in / 50.8 cm	18.5 in / 46.9 cm	32 in / 81.2 cm	38 in / 96.5 cm	300 lb / 136 kg
716170A-4	17 in / 43.1 cm	22 in / 55.8 cm	20 in / 50.8 cm	18.5 in / 46.9 cm	25 in / 63.5 cm	32 in / 81.2 cm	300 lb / 136 kg
716270A-1 716270A-2 716270A-4 716270 ColorSelect series	17 in / 43.1 cm	25.5 in / 64.7 cm	20 in / 50.8 cm	20.2 in / 51.3 cm	32 in / 81.2 cm	39 in / 99 cm	300 lb / 136 kg
716370A-2 716370A-4	17 in / 43.1 cm	24 in / 60.9 cm	20 in / 50.8 cm	18 in / 45.7 cm	25.5 in / 64.7 cm	32.5 in / 82.5 cm	300 lb / 136 kg
613070A	17.25 in / 43.8 cm	23.2 in / 58.9 cm	19.7 in / 50 cm	19.5 in / 49.5 cm	32.5 in / 82.5 cm	39.5 in / 100.3 cm	300 lb / 136 kg
700175C-2 700175CR	18 in / 45.7 cm	26 in / 66 cm	23 in / 55.4 cm	20 in / 50.8 cm	32 in / 81.2 cm	39 in / 99 cm	400 lb / 181.4 kg
1925A	17 in / 43.1 cm	22.5 in / 57.1 cm	19.7 in / 50 cm	20.5 in / 52 cm	31 in / 78.7 cm	38 in / 96.5 cm	400 lb / 181.4 kg
6015A	-	16.5 in / 41.9 cm	-	13.2 in / 33.5 cm	29.5 in / 74.9 cm	34.5 in / 87.6 cm	300 lb / 136 kg

ROLLATORS

The Lumex rollators are designed to allow for easy maneuverability and comfort to aid in mobility. Each rollator frame is built for durability and features the following:

- Padded seat for added comfort
- Easy to remove padded backrest in standard and contoured style
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Adjustable handle height to accommodate a broad range of user heights
- Comes with wire basket or pouch for storage of personal items



E0143*
E0156

Walkabout Lite Four-Wheel Rollator

- Lightweight aluminum frame weighs only 14 ½ lb
- 6" wheels designed for easy turning and mobility
- Basket for storage can be placed under seat or in front of frame
- Reinforced frame built for durability
- Available in 8 attractive, retail-friendly color options
- Item No. RJ4300B, RJ4300G, RJ4300K, RJ4300L, RJ4300P, RJ4300R, RJ4300AQ, RJ4300CH

Walkabout Wide Four-Wheel Rollator

- Wider, aluminum frame with 18.5 inch width in between handgrips allows for additional patient comfort when positioning and seated for resting
- Wide 15.5 inch (39.3 cm), padded seat provides larger seating surface for comfort
- 6 inch wheels designed for easy turning and mobility
- Item No. RJ4318B, RJ4318K, RJ4318R, RJ4318AQ



E0143*
E0156

Walkabout Basic Four-Wheel Rollator

- Lightweight 13 ½" aluminum frame is an affordable solution to mobility yet maintains ease of use for patient
- Convenient black storage pouch for personal items
- 6" wheels designed for easy turning and mobility
- Item No. RJ4900B, RJ4900R



E0143*
E0156

Walkabout Essentials Four-Wheel Rollator

- Durable, cost effective steel frame
- Padded seat for added comfort
- Removable, hinged, padded backrest can be folded for transport
- 6" wheels designed for easy turning and mobility
- Pouch included for storage of personal items
- Item No. RJ5500R, RJ5500K



E0143*
E0156

Walkabout Junior Four-Wheel Rollator

- Smaller frame is ideal for petite or junior users
- 6" wheels designed for easy turning and mobility
- Tool-free removable padded backrest for comfort
- Item No. RJ4301R



E0143*
E0156

*HCPCS Code



colorSelect
Available on select models



E0143*
E0156

Walkabout Deluxe Contour Four-Wheel Rollator

- Contoured, padded backrest offers added comfort
- 8" wheels are ideal for indoor and outdoor use
- Comes with convenient basket placed under seat or in front frame for personal items
- Reinforced frame is built for durability and long term use
- Padded seat gives user comfort when seated
- Item No. RJ4805B, RJ4805K, RJ4805L, RJ4805P, RJ4805R, RJ4805CH

Walkabout Hemi Four-Wheel Rollator

- Lower seat-to-floor height accommodates shorter individuals
- Convenient wire basket for personal items
- 6" wheels designed for easy turning and mobility
- Item No. RJ4302B, RJ4302K, RJ4302R



E0143*
E0156



E0143*
E0156

Set N' Go Height Adjustable Rollator

- Offers the convenience and flexibility of a two-in-one product all rolled into one
- Easy seat height adjustment to be positioned from 18" to 22" in 1" increments
- Comfortable, seamless padded seat with zippered pouch under seat and larger black pouch
- Adjustable handle height with knobs to secure in place
- 6" wheels designed for easy turning and mobility
- Item No. RJ4700B, RJ4700R, RJ4700S, RJ4700P

Set N' Go Wide Height Adjustable Rollator

- Offers the same features as the standard Set N' Go but in a Wide option with 19" width in between handgrips
- Easy height adjustment on the frame allows seat height to be positioned from 19.5" to 23.5" in 1" increments
- Wider 15.75" comfortable, seamless padded seat with zippered pouch under seat for security of personal items
- Also features Velcro® folding strap to keep rollator folded for storage and transport
- Item No. RJ4718B, RJ4718R, RJ4718S



E0143*
E0156

Walkabout Imperial Bariatric Four-Wheel Rollator

- Black 8" wheels are ideal for indoor and outdoor use
- Padded oversized 20" seat
- Wide frame made of high quality aluminum to remain lightweight
- Comes with large pouch under the seat for storage of personal items
- Easy to remove padded backrest for comfort
- Item No. RJ4400R, RJ4400K (Standard Backbar)
- Item No. RJ4405B, RJ4405K, RJ4405R (Contoured Backbar)

Walkabout Imperial Hemi Bariatric Four-Wheel Rollator

- Same features as the standard Bariatric but in a hemi 19" seat-to-floor height for shorter individuals
- Item No. RJ4402R



E0149*
E0156



Three-Wheel Cruiser

- Provides convenient mobility for those not needing the extra stability of a 4 wheel rollator
- Lightweight aluminum frame
- Available in 3 attractive colors: Pewter Metallic, Metallic Burgundy and Bondi Blue
- Adjustable handle height for broad range of users
- 8" black wheels for indoor and outdoor use
- Comes complete with zippered carry pouch
- Item No. 609201B, 609201M, 609201P

E0143*

HybridLX Rollator/Transport Chair

- Lightweight, durable aluminum frame weighs only 22 lbs (this includes footrests)
- Hybrid concept – designed to be used primarily as a rollator while easily transforming into a transport chair with removable footrests (adjustable in length) when user needs a short rest
- Ultra-Compact design – easy, tool-free assembly and folds effortlessly for convenient transport and storage
- Flexible, soft polyester backrest and 18" wide mesh seat with pull up strap for easy folding provides comfort for user
- Large 8" wheels for easy rolling on smooth and other terrain
- Item No. LX1000T, LX1000B



E0143*
E0156



SureStride Rollator

- Patented seat slides along frame allowing user to walk inside the rollator frame allowing for improved gait
- Adjustable push bar handle provides a place to propel rollator forward with better posture and stability
- Ergonomic handgrips with easy brake access
- Intuitive one brake handle can activate both brakes
- Deluxe, soft pouch for personal items can be easily installed and removed
- Large wheels allow for easy maneuverability
- Item No. S9000TS

Item No.	Width Inside Handgrips	Overall Width	Overall Depth	Seat Depth	Seat Width	Seat Height	Handle Height Min/Max	Weight Capacity ¹ EVENLY DISTRIBUTED
RJ4300B, G,K,L,P,R,AQ,CH	16 in / 40.6 cm	24.5 in / 60.9 cm	27.5 in / 69.8 cm	14.5 in / 36.8 cm	13.5 in / 34.2 cm	23 in / 55.4 cm	32 in - 37 in 81.2 cm - 93.9 cm	300 lb / 136 kg
RJ4318B, G,K,L,P,R,AQ,CH	18.5 in / 46.9 cm	27.5 in / 69.8 cm	27.5 in / 69.8 cm	15.5 in / 39.3 cm	13 in / 33 cm	22.5 in / 57.1 cm	32 in - 37 in 81.2 cm - 93.9 cm	350 lb / 158.7 kg
RJ4900B,R	16.5 in / 41.9 cm	24.2 in / 61.4 cm	24.5 in / 62.2 cm	13 in / 33 cm	12.5 in / 31.7 cm	21 in / 53.3 cm	32.5 in - 37.5 in 82.5 cm - 95.2 cm	250 lb / 113.3 kg
RJ4805B, G,K,L,P,R,CH	16 in / 40.6 cm	25 in / 63.5 cm	27.5 in / 69.8 cm	14 in / 35.5 cm	13.5 in / 34.2 cm	22.5 in / 57.1 cm	34 in - 38.5 in 86.3 cm - 97.7 cm	300 lb / 136 kg
RJ5000	15.5 in / 39.3 cm	24 in / 60.9 cm	30.7 in / 77.9 cm	12 in / 33.4 cm	13 in / 33 cm	19.7 in / 50 cm	30 in - 36 in 76.2 cm - 91.4 cm	300 lb / 136 kg
RJ4301R	12 in / 33.4 cm	20 in / 50.8 cm	24 in / 60.9 cm	12 in / 33.4 cm	12 in / 33.4 cm	18.5 in / 46.9 cm	27.5 in - 32 in 69.8 cm - 81.2 cm	300 lb / 136 kg
RJ4302B,K,R	15.5 in / 39.3 cm	22.5 in / 57.1 cm	21.5 in / 54.6 cm	14 in / 35.5 cm	13.5 in / 34.2 cm	19.5 in / 49.5 cm	32 in - 36.5 in 81.2 cm - 92.7 cm	300 lb / 136 kg
RJ4700B,R,S	16.3 in / 41.4 cm	24.5 in / 60.9 cm	29 in / 73.6 cm	14 in / 35.5 cm	12 in / 33.4 cm	18" - 22"	29 in - 38 in 73.6 cm - 96.5 cm	300 lb / 136 kg
RJ4718B,R,S	19.4 in / 49.2 cm	26.9 in / 68.3 cm	27.5 in / 31.5 in 69.8 cm / 80 cm	12 in / 33.4 cm	15.7 in / 39.8 cm	19.5" - 23.5"	32 in - 38 in 81.2 cm - 96.5 cm	350 lb / 158.7 kg
RJ4400K,R, RJ4405B,K,R	24 in / 60.9 cm	29.5 in / 74.9 cm	27.5 in / 69.8 cm	20 in / 50.8 cm	13.5 in / 34.2 cm	23.5 in / 59.6 cm	34 in - 39 in 86.3 cm - 99 cm	450 lb / 204.1 kg
RJ4402R	22.5 in / 57.1 cm	29.5 in / 74.9 cm	30 in / 76.2 cm	20 in / 50.8 cm	13.5 in / 34.2 cm	20.5 in / 52 cm	31 in - 34 in 78.7 cm - 86.3 cm	450 lb / 204.1 kg
LX1000B,T	18.5 in / 46.9 cm	26 in / 66 cm	40 in / 101.6 cm	28.5 in / 52.1 cm	18 in / 45.7 cm	300 lb / 136 kg	30 in - 37 in 76.2 cm - 93.9 cm	300 lb / 136 kg
S9000TS	16 in / 40.6 cm	24.5 in / 60.9 cm	31.7 in / 80.5 cm	15.1 in / 38.4 cm	17.2 in / 43.7 cm	22 in / 55.9 cm	29.4 in - 35.4 in 74.7 cm - 89.9 cm	300 lb / 136 kg
609201B,M,P	14.5 in / 36.8 cm	25 in / 63.5 cm	22.5 in / 57.1 cm	-	-	-	33 in - 37 in 83.8 cm - 93.9 cm	250 lb / 113.3 kg

Color Select is a registered trademark of the Rollex Corporation

¹ maximum weight capacity - evenly distributed - refer to user manual *HCPCS Code

GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

published **June 2026**

Images are digital reproductions. Perceived texture and coloring may vary on
finished product. Monitor and printer settings will effect appearance of color.

+1 770.368.4700

Information contained herein is subject to change without notice. The most
current & complete product information can be found on our website.

www.grahamfield.com



© 2014, GF Health Products, Inc. All Rights Reserved
Lumex, is a registered trademark of GF Health Products, Inc.
Color Select is a registered trademark of the Rollex Corporation

graham
field

GF400050RevF26