

## MEDICAL SEATING

# PHYSICIAN EXAM STOOL

## 9000 SERIES

Lumex® offers a variety of medical stools to meet the needs of an active facility. Manufactured with the highest quality materials to ensure a long life, Lumex stools offer the most wanted features in an exam stool.

### Features:

- Seat diameter comfortably fits to most users
- Pneumatic assist height adjustment for ease of use
- Choice of 3 base styles to complement facility decor
- Choice of height adjustment lever options
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate
- Wide diameter base with hard wheel casters for greater stability
- Backrest, seat height, and footring options
- Available in 24 industry standard colors to complement most decors

Product Info 9501 9702 9903 Fabric Choices



### 9501 Options:

black nylon composite base +

|               | Standard Height | with Backrest | with Footring | Lab Height |
|---------------|-----------------|---------------|---------------|------------|
| 9501-xx       | ■               |               |               |            |
| 9501-B-xx     | ■               | ■             |               |            |
| 9501-FR-xx    | ■               |               | ■             |            |
| 9501-BFR-xx   | ■               | ■             | ■             |            |
| 9501-10BFR-xx |                 | ■             | ■             | ■          |

★ options available at additional cost



### 9702 Options:

brushed aluminum base and toe caps +

|               | Standard Height | with Backrest | with Footring | Lab Height |
|---------------|-----------------|---------------|---------------|------------|
| 9702-xx       | ■               |               |               |            |
| 9702-B-xx     | ■               | ■             |               |            |
| 9702-FR-xx    | ■               |               | ■             |            |
| 9702-BFR-xx   | ■               | ■             | ■             |            |
| 9702-10BFR-xx |                 | ■             | ■             | ■          |

★ options available at additional cost



### 9903 Options:

polished aluminum base +

|              | Standard Height | with Footring | Lab Height |
|--------------|-----------------|---------------|------------|
| 9903-xx      | ■               |               |            |
| 9903-FR-xx   | ■               | ■             |            |
| 9903-10FR-xx |                 | ■             | ■          |

★ options available at additional cost

### Specifications:

| Series | Seat Diameter |       | Base Diameter with casters |       | Standard Seat Height |                    |                    |                    | Lab Seat Height |       |       |       | Overall Height w/ Backrest |                    |                    |                    | Footring Diameter |       | Caster Diameter hard wheel |      | Shipping Weight <sup>8</sup> |                 | Weight Capacity <sup>1</sup> EVENLY DISTRIBUTED |                 | Pneumatic Height Adjustment | Single Lever D-Handle | 360° Chrome Handle |
|--------|---------------|-------|----------------------------|-------|----------------------|--------------------|--------------------|--------------------|-----------------|-------|-------|-------|----------------------------|--------------------|--------------------|--------------------|-------------------|-------|----------------------------|------|------------------------------|-----------------|-------------------------------------------------|-----------------|-----------------------------|-----------------------|--------------------|
|        | inch          | cm    | inch                       | cm    | low                  | high               | low                | high               | low             | high  | low   | high  | low                        | high               | low                | high               | inch              | cm    | inch                       | cm   | lb                           | kg <sup>1</sup> | lb                                              | kg <sup>1</sup> |                             |                       |                    |
| 9501   |               |       |                            |       | 18.50                | 46.99              | 24.00              | 60.96              | 24.50           | 62.23 | 35.00 | 88.90 | 33.00                      | 83.82              | 39.00              | 99.06              |                   |       |                            |      | 20                           | 9.07            |                                                 |                 | ■                           | ■                     |                    |
| 9702   | 16.00         | 40.64 | 24.00                      | 60.96 | 17.50                | 44.45              | 23.00              | 58.42              | 23.50           | 59.69 | 34.00 | 86.36 | 32.50                      | 82.55              | 37.50              | 95.25              | 18.00             | 45.72 | 2.00                       | 5.08 | 22                           | 9.98            | 350                                             | 158.76          | ■                           |                       | ■                  |
| 9903   |               |       |                            |       | 18.00 <sup>2</sup>   | 45.72 <sup>2</sup> | 23.00 <sup>2</sup> | 58.42 <sup>2</sup> | 22.00           | 55.88 | 32.00 | 81.28 | 18.00 <sup>2</sup>         | 45.72 <sup>2</sup> | 23.00 <sup>2</sup> | 58.42 <sup>2</sup> |                   |       |                            |      | 23                           | 10.43           |                                                 |                 | ■                           |                       | ■                  |

★ options available at additional cost † -xx in the product number is the placeholder for fabric choice ‡ standard rounding method used  
<sup>1</sup> maximum weight capacity - evenly distributed - refer to user manual <sup>2</sup> overall height without backrest <sup>8</sup> approximate measurement



GF Health Products, Inc. is an  
 intertek ISO 13485:2016 / MDSAP Certified Company

•1 770.368.4700

Information contained herein is subject to change without notice. The most current & complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022 GF Health Products, Inc. All Rights Reserved  
 Lumex is a registered trademark of GF Health Products, Inc.  
 Lumex, Everest & Jennings, and Graham Field are trademarks of GF Health Products, Inc.

**graham field**