

Universal Raised Toilet Seat with Armrests

- Built in toilet safety rail.
- Textured arm rest designed to provide good gripping surface.
- Seat height adjusts from 3 1/2" to 6 1/2".
- Arm height above seat: 7 1/2".
- Full length splash guard extends into toilet for greater protection.
- Easily removed when not in use.
- White plastic seat is chip and crack resistant.
- Fits most toilets.
- 300 lb maximum weight capacity.
- Designed to meet the requirements of HCPCS code: E0244.
- Lifetime limited warranty.

6990A 1 ea



ITEM	Width Between Arms	Overall Width	Seat Width	Seat Depth	Weight Capacity	Replacement Tip
6990A	17 1/2"	19 1/2"	14"	16"	300 lb	NA



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com